

System requirements:

- for Matlab Code:
 - Windows 7, 64 bit or later
 - Matlab v. 2017b or later with "Image Processing and Statistics" and "Machine Learning" toolbox
 - No specific hardware required
- for compiled Software:
 - Windows 7, 64 bit or later
 - Matlab Runtime v. 2018a
 - No specific hardware required

Installation guide:

- Copy code to Matlab Path, download dependencies (see Wiki on <https://github.com/DWALab/Phindr3D>) and run from Matlab
- or**
- Install from the provided EXE file and install Matlab runtime if required.

Demo:

To test the functionality of Phindr3D please download an excerpt of the full data sets from Zenodo: The sample data include:

- neuron data set to test the functionality of all Phindr3D core functions (20 3-channel image stacks from four experimental groups from one experimental replicate; size of the zipped data set: 13.8 GB; <https://dx.doi.org/10.5281/zenodo.4064148>)
- organoid data set to test the functionality of the organoid selection and cropping feature (4 3-channel image stacks with 81 z-planes; size of the zipped data set: 5.4 GB; <https://dx.doi.org/10.5281/zenodo.4384912>).

Expected output:

- Neuron data set: visualization of the feature space the can be changed using PCA, t-SNE or Sammon mapping; further analysis of the data using clustering or classification. All output plots can be saved from the Phindr3D GUI.
- Organoid data set: cropped organoid volumes will be saved as individual TIF files into the specified output folder.


Expected run time on a "normal" desktop computer:

- neuron data set: ~30 minutes to retrieve feature data
- organoid data set: ~20 minutes to crop all organoids

Instructions for Use


MAIN GUI

- Go to folder where code resides. Type <Phindr3D> from Matlab command prompt
- The following window opens up:



CREATE METADATA


- Click on Metadata from Menu > Metadata > Create Metadata



- Select root directory where images are located (Red box)
- The sample file name automatically fills up
- Type regular expression to extract metadata from file name (Green box).
- *When using the test data or other data generated using PerkinElmer Harmony, type:*
`(?<Well>\w+)f(?<Field>\d+)p(?<Stack>\d+)-ch(?<Channel>\d).*.tif(f)?`

- For a tutorial see https://www.mathworks.com/help/matlab/matlab_prog/regular-expressions.html
- You can now evaluate the regular expression
- Write a prefix for the file name (Yellow box) and press “Create File”. A metadata file is created in the same directory as the images
- Close this window manually


LOAD METADATA



- Press Load Metadata and select the metadata file.
- Select channels in the new window and press ok

SET PARAMETERS

- Standard parameters are set automatically and work fine for many situations.
- However, if the user needs to change parameters press “Set Voxel Parameters”. A new window pops up where the parameters can be changed.




The "Set Parameters" dialog box is a window for configuring analysis parameters. It features two main sections: "Super Voxel" and "Mega Voxel". Each section contains a "Size" sub-section with input fields for X, Y, and Z dimensions. Below these are input fields for the number of categories: "#SV Categories" (set to 15) and "#MV Categories" (set to 40). Further down, there are fields for "#Voxel Categories" (set to 20) and "#Training Images" (set to 5). There are also checkboxes for "Use background voxels for computation", "Normalise Intensity Per Condition", and "Train by condition". At the bottom, there are two "Plate" dropdown menus and "Reset" and "Done" buttons.

Parameter	Value
Super Voxel Size X	10
Super Voxel Size Y	10
Super Voxel Size Z	3
#SV Categories	15
Mega Voxel Size X	5
Mega Voxel Size Y	5
Mega Voxel Size Z	2
#MV Categories	40
#Voxel Categories	20
#Training Images	5


- Press “DONE”

RUNNING PHINDR3D

- Multi-channel image is displayed in the main window
- Press “Next Image” to view another image in the folder




- To change channel colors Press Channel colors



- Click on individual channels and select the color for each channel. Close the window once you set colors
- To see SV/MV, check the SV or MV box
- To extract features, Press **Phind** button. Enter file name for the output of the features
- Once completed, a feature file is created in the output directory provided by the user

VIEW PHINDR3D RESULTS

- Click <View Results> from the main window




- Select output feature file by clicking either <Select Feature File> or File>> Input Feature File from the main menu (Red).

- Once the feature file is selected, the software will ask the user to define Metadata (i.e. all columns except the features [MV1, MV2, etc]) from the feature file. After this, the software will ask to select the Channel columns.
- Phindr3D performs a PCA and loads the plot on the plot window. The user can change the type of plot using the drop down menu as well as choose between 2D/3D plots (Brown box)
- The user can color the points using one of the metadata Columns (Green box)
- The plots can be exported as vector graphics (SVG file)

SHOW IMAGE

- To see the corresponding image for any point on the plot, hover over the point, right click and click show image. A new image window is shown:



CLUSTER IMAGE DATA

- *Automated Cluster Estimation*: Click Data Analysis >> Clustering >> Estimate Clusters
- *Manual Clustering*: Click Data Analysis >> Clustering >> Set Number of Clusters
- *Pie Charts*: Click Data Analysis >> Clustering >> Pie Maps (Only available after clustering)
- *Export*: Click Data Analysis >> Clustering >> Export Clustering Results (Only available after clustering)

CLASSIFY IMAGE DATA

- Click Data Analysis >> Classification >> Select Classes
- Select the classes and click OK. Once the classifier is built, the output is saved at a user-defined location